


Practical
and well prepared



VERTEBRIS Sieve Basket Systems

Sterilization and storage
of the **VERTEBRIS** Instrument Set


VERTEBRIS Sieve Basket Systems

Sterilization and storage in conformity with the latest standards

VERTEBRIS sterilization baskets with logical arrangement of the instruments permit optimized reprocessing of your valuable **VERTEBRIS** Instrument Set and always ensure a better overview and controlled access.

The various levels are printed with the contours of the individual instruments. This simplifies loading while also ensuring scope for individualized instrument configuration.

- Proven Richard Wolf Sterilization Sieve Basket System
- Optimum sterilization with individualized point support
- Stable design with patented double frame



Sieve Basket System for the VERTEBRIS lumbar-thoracic Instrument Set



VERTEBRIS lumbar-thoracic Instrument Sieve Basket
for accommodating
endoscopes, access and working instruments,
accessories for radio-frequency surgical system
L x W x H: 530 x 250 x 150 mm 85841217

Sieve Basket System for the VERTEBRIS stenosis Instrument Set



VERTEBRIS stenosis Instrument Sieve Basket
for accommodating
endoscopes, access and working instruments,
accessories for radiofrequency surgery and motor system
L x W x H: 530 x 250 x 150 mm 85841223

The instrument sieves shown above are supplied without instruments.
You can obtain detailed information on **VERTEBRIS** Instrument Sets for full-endoscopic spine surgery by contacting the Sales Field Service.