



Aequalis™ Cemented

Your experience guides our choices

Developed in collaboration with the greatest names in shoulder surgery, Aequalis™ is the first third-generation modular and adaptable shoulder prosthesis.

CEMENTED*

Adaptability for reproducing movement

HUMERAL COMPONENT

Inclination of the humeral head

- Reproduces the inclination of the patient's humeral head
- 4 angles of inclination: 125°, 130°, 135°, 140°

Stem design

- Optimized positioning
- 4 diameters: 6.5 mm, 9 mm, 12 mm, 15 mm

Variable retroversion

- Simple and surgeon-friendly instrumentation helps to reproduce intra-operatively each patient's exact anatomical retroversion.

Offcenter dial

- Enables the humeral head to be offset to reproduce each patient's anatomical characteristics.

Head design

- Link between diameter and thickness of the humeral head for optimum and anatomical joint restoration.

Designed to reproduce anatomic characteristics

- Reconstructs joint surfaces
- Restores anatomical center of rotation
- Respects joint biomechanics

* Maitrise Orthopédique #26 - September 1993. Adaptabilité et modularité au cours des prothèses d'épaule. P. Boileau, Centre Hospitalier Universitaire de Nice 06006 Nice Cedex1 ; G. Walch, Centre Hospitalier Universitaire de Lyon 69310 Pierre-Bénite.



IMPLANTS

Cemented humeral stem

	Ø	125°	130°	135°	140°
Standard	6.5 mm	DWB061	DWB062	DWB063	DWB064
	9 mm	DWB091	DWB092	DWB093	DWB094
	12 mm	DWB121	DWB122	DWB123	DWB124
	15 mm*	DWB151	DWB152	DWB153	DWB154
Long	6.5 mm*		DWB562		
	9 mm*		DWB592		
	12 mm*		DWB622		
Extralong	9 mm*		DWB662		
	12 mm*		DWB663		
	15 mm*		DWB664		

*upon request only

Humeral head

Size	Ref
37 x 13.5 mm*	DWB237
39 x 14 mm	DWB239
41 x 15 mm	DWB241
43 x 16 mm	DWB243
46 x 17 mm	DWB246
48 x 18 mm	DWB248
50 x 16 mm	DWB250
50x 19 mm	DWB251
52x 19 mm*	DWB252
52 x 23 mm*	DWB253
54 x 23 mm*	DWB254
54 x 27 mm*	DWB255

*upon request only

Keeled glenoid

Description	Ref
S30 Keeled Glenoid	DWE501
S35 Keeled Glenoid	DWE502
S40 Keeled Glenoid	DWE503
M30 Keeled Glenoid	DWE511
M35 Keeled Glenoid	DWE512
M40 Keeled Glenoid	DWE513
L40 Keeled Glenoid	DWE521
L50 Keeled Glenoid	DWE522
L60 Keeled Glenoid	DWE523
XL40 Keeled Glenoid	DWE531
XL50 Keeled Glenoid	DWE532
XL60 Keeled Glenoid	DWE533

Cement restrictor

Size	Ref
Humerus / Ø 24 mm	EBO101
Ulna / Ø 13 mm	EBO102


US HEADQUARTERS

Tornier, Inc.
10801 Nesbitt Avenue South
Bloomington, MN 55437
USA
+1 952 426 7600

INTERNATIONAL HEADQUARTERS MANUFACTURER

Tornier SAS
161 rue Lavoisier
38330 Montbonnot Saint Martin
FRANCE
+33 (0)4 76 61 35 00

www.tornier.com

Prior to using any Tornier device, please review the instructions for use and surgical technique for a complete listing of indications, contraindications, warnings, precautions, potential adverse events, and directions for use. © 2014 Tornier, Inc. All rights reserved. Aequalis™, Tornier™ and  are trademarks or registered trademarks of Tornier in the U. S. and other countries.

Aequalis™ Cemented - Brochure - UDAF151