


PHALINX[®]

Hammertoe Fixation



PHALINX[®]



Hammertoe Fixation

System Advantages

- Titanium implants eliminate need for freezer storage
- Implants available in straight cannulated and 10° angled
- Color-coded implants and instruments increase system efficiency



Simple & Effective Hammertoe Fixation



| Implant Description / Size | Length (mm) | | | Diameter (mm) | | |
|-------------------------------|-------------|------------|----------|---------------|----------|-------------|
| | A Overall | B Proximal | C Distal | D Proximal | E Distal | F Thickness |
| Extra Small Cannulated/Angled | 11.2 | 6 | 4.2 | 2.5 | 3.2 | 1.8 |
| Small Cannulated/Angled | 13.0 | 7 | 5.5 | 2.8 | 4.0 | 1.9 |
| Medium Cannulated/Angled | 14.1 | 7 | 6.0 | 3.0 | 4.5 | 2.0 |
| Large Cannulated/Angled | 15.2 | 8 | 6.2 | 3.3 | 5.0 | 2.0 |



Wright Medical Technology, Inc.
1023 Cherry Road
Memphis, TN 38117
800 238 7117
901 867 9971
www.wmt.com

Wright Medical EMEA
Atlas Arena, Australia Building
Hoogoorddreef 7
1101 BA Amsterdam
the Netherlands
011 31 20 565 9060