

ORTHOLOC® 3Di

Foot Reconstruction System

BOW® Plate

Anatomic 1st metatarsal opening wedge plate for reliable hallux valgus correction



- **Anatomic Design:** Low profile anatomic plate contours designed to reduced soft tissue irritation
- **Versatility:** 2.7 or 3.5mm polyaxial locking screws provide reliable construct strength
- **Controlled Correction:** 0-5mm wedge spacers allow for precise IM angle correction

More Information:

- ORTHOLOC® 3Di Foot Reconstruction System Hallux Surgical Technique FA285-412
- ORTHOLOC® 3Di Foot Reconstruction System BOW® Plate Animation FA290-412

Additional Products:

- DARTFIRE® Small Screw System
- CHARLOTTE® Quick Staple
- ORTHOLOC® 3Di 1st Metatarsal Plate



 **WRIGHT**[™]
FOCUSED EXCELLENCE