


# ORTHOLOC® 3Di

## Ankle Fusion Plating System

### Options for Each Patient



- **Patient-Specific Fixation:** plate holes accept both 4.5mm and 5.5mm non-locking and variable angle locking screws
- **Controlled Compression:** plate compression slots allow for eccentric screw placement and compression, while the ORTHOLOC® 3Di X-Track provides additional compression external to the plate
- **Comprehensive Implant Offering:** system features 6 different anatomical plates designed to facilitate either an anterior, lateral or posterior approach

#### More Information:

- ❑ ORTHOLOC® 3Di Ankle Fusion Plating System Surgical Technique 007794
- ❑ ORTHOLOC® 3Di Ankle Fusion Plating System Brochure 007791
- ❑ ORTHOLOC® 3Di Ankle Fusion Plating System Animation 007386

#### Additional Products:

- ❑ VALOR® Hindfoot Fusion Nail System
- ❑ DARCO® Headed Screws

